

HPLC/GC Accessories Book



LC Columns

Luna C8, C18, Silica(2), CN, NH2, SCX, HILIC, PFP(2)

Kinetex C8, C18, XBC18, BiphenylPhenylHexyl, F5

Gemini C6, C18, NXC18, C6Phenyl

Jupiter C4, C5, C18, Proteo

Rezex RHM, RCM, ROA, RSO, RNO, RAM

Lux Chiral Cellulose1, Cellulose2, Cellulose3,

Cellulose4, Amylose1, Amylose2, Cellulose5

Synergi MAXRP, HydroRP, Polar RP, Fusion RP

GC Columns (Zebron)

ZB1, ZB1ms,

ZB1HT Inferno, ZB1XT SimDist,

ZB5, Zb5Msi, ZB5HT Inferno,

Zb5Ms, ZB 35, ZB5MSplus,

ZB SemiVolatiles, ZB35HT

Inferno, ZB50, ZB624, ZB1701,

ZB1701P, ZBWAX, ZBWAXplus



Zebron



Gemini

SecurityGuard



KINETEX




LUX



LUNA

ONYX



REZEX



PHENEX
Syringe Filter



EZfaast™
Amino Acid Analysis Kits



strata-X

Solid Phase
Extraction (SPE)



ชุดวิเคราะห์กรดอะมิโน



roQ QuEChERS Kit



VEREX
Certified Vial Products
Premium Grade
Vials



ชุดเตรียมตัวอย่าง

LC Columns

Hypersil BDS C8, C18, Phenyl, ODS, MOS(C8),
Hypersil Gold C8, C4, aQ, PFP, Phenyl, CN, Silica, HILIC
Accliam 120 C8, 120 C18, Phenyl 1, C30, Amg C18,
 Polar Advantage, Polar Advantage II
Accucore RPMS, C8, C18, aQ, Phenyl – Hexyl, PFP, Polar
 Premium, C30, HILIC, XL C18, XL C8
Synchronis C8, C18, aQ, Phenyl, Amino, Silica, HILIC
Accucore Vanquish

GC Columns

Trace TR1, TR5, TRWAX, TR35MS, TRFFAP
Trace Gold TG1MS, TG5MS, TG5HT,
 TGSQC, TG17MS, TG225MS, TGVMS
Trace Plot TGBOND Q, TGBOND Alumina,
 TGBOND Msieve 5A, TGBOND U



Trace™



Hypersil BDS™



Hyper Carb™



Hypersil Gold™



Synchronis™



Accucore™



Syringe Filter



Tip Microscale SPE
Extraction Tips, 110µL



Solid Phase
Extraction (SPE)



QuEChERS Kit

ชุดเตรียมตัวอย่าง



Premium Grade
Vials

LC Columns

ACE Excel NH2, AQ, C18, C18-AMIDE, C18-AR, C18-HL, C18-PFP, C4, C8, CN, CN-ES, PHENYL, SIL, SUPER C18

ACE UltraCore SUPERC18, SUPERPHENYLHEXYL


ACE AQ, C18, C18-300, C18-AMIDE, C18-AR, C18-PFP, C18-HL, CN-ES, NH2, C18HL, C8, C8-300, C4, C4-300, CN, CN-300, PHENYL, AQ, SIL, SUPERC18



ACE Super C18

ACE Excel

ACE Ultra Core



Coverage for Various kinds of Applications



LC Columns

Alltima Amino, C18, C18-LL, C8, Cyano, Phenyl, Silica

Alltima HP C18, C18-Amide, C18-AQ, C18-EPS, C18-HL, C8, CN, HILIC, Silica

Apollo C18, C8, Phenyl, Silica

Previal Amide, Amino, C8, C18, Carbohydrate ES, Cyano, Organic Acid, Silica

Apex II Amino, II C18, ODS, Phenyl, Phenly-RP, Silica

Allsep Anion IC

Genesis C18, AQ, C8, Cyano, Phenly, Silica, C4-300


Vydac 214MS, 214TP, 218MS, 218TP, DENALI, EVEREST

GRACE

USA

www.grace.com





LC Columns

Pack Pro C8, C18

Pack ODSA, ODSAM, ODSAQ, ODSAL, C4, Ph, CN, Hydrosphere

CHIRAL CD BR, PREP CD ST/PM, NEA(R) (S)

CHIRAL ART AmyloseC, CelluloseC, AmyloseSA, CelluloseSC, AmyloseSE, Alcyon SFC

Triart C8, C18, C18 ExRS, Phenyl

Carotenoid

J'sphere ODSH80, ODSM80, ODSL80, C18

UltraHT Pro C18, Hydrosphere C18

Meteoric Core C18, C18 BIO

LC Columns

Europa C 4, C8, C18

Brisa LC2 C18

Tracer Excel 120 Silica, C1, C4, C8, ODS-A, ODS, APS, CN, P

Tracer Extrasil C-1, C-8, C-6, ODS1

Mediterranea Sea 4, Sea 8, Sea 18

Hyperpack BASIC, C8, ODS



VDSpher – the material!



High-End End-Capping


- excellent batch-to-batch reproducibility
- large selection of RP phases
- large selection of particle sizes (1.8 µm to 30 µm)
- very high degree of purity



Quality Made in Germany since 1987



Compounds (1) Phenol (2) Pyridine



Columns

VDSpher® (Classic) 100 C18-E
 VDSpher® PUR 100 C18-E
 (5 µm, 250 mm x 4.6 mm)

Chromatographic conditions

Solvents: A: Acetonitrile
 B: Water
 Method: isocratic
 50 % A/50 % B
 Flow rate: 1.0 ml/min
 Pressure: VDSpher® (Classic): 10.3 MPa
 VDSpher® PUR: 9.6 MPa
 Temp.: 30 °C
 Detection: UV (λ = 254 nm)

LC Analyticals Columns

VDSpher (Classic)

SIL, C4, C8, C18, Phenyl, CN, Diol, NH₂

VDSpher PUR

SIL, C4, C8, C18, Phenyl, CN, Diol, SAC, NH₂

U-VDSpher PUR SIL, C8, C18, CN, Phenyl

VDSpher Opti Aqua C4, C8, C18

VDSpher Opti Bio C4,C8, C18

Size Specification

Porous Size (Å) : 75, 100, 150, 200, 300
 Particle Size (µm): 1.8, 2.2, 3, 3.5, 4, 5, 10
 Diameter(mm): 1, 2, 3, 4, 4.6
 Length (mm): 30, 33, 50, 53, 75, 100, 125,
 150, 200, 250, 300

LC Semi-Preparative Columns

VDSpher (Classic)

SIL, C4, C8, C18, Phenyl, CN, NH₂

VDSpher PUR SIL, C4, C8, C18

VDSpher Opti Aqua C18

Size Specification

Porous Size (Å) : 75,100
 Particle Size (µm) : 5,10 µm
 Diameter(mm) : 8,10,16,20 mm
 Length (mm) : 30,40,50,60,100,125,150,250

Other Accessories

Syringe Filter
 Membrane
 Filtration Apparatus
 Hamilton Syringe
 Rubber Free Syringe
 Vial
 Tape
 Microwave Vessel

HAMILTON
 THE MEASURE OF EXCELLENCE™



Coloured Labelling Tape



Used in the 121 °C moist heat sterilization



Microwave digestion vessel for milestone & CEM



Water purification systems/VWR
 From tap Type 2 or Type 1 ultrapure water
 Produce up to 60 L/hr.



CombiFlash Rf+
 Sample 4-330 g.
 Flowrate 1-200 ml/min



Membrane Solutions

Lubitech Technologies



www.membrane-solutions.com www.lubitech.com

www.anpel.com.cn

<http://jinteng1993.com>

Our Products

1. Young Lin Instrument (Korea) - HPLC, GC, GC-MS, LC-MS, Water Purification System
2. Lab Tech (USA.) - Sample Preparation System
3. Data Apex (The Czech Republic) - Clarity Chromatography Software
4. Pickering Laboratories (U.S.A.) - Post-Column Derivatization, Amino Acid
5. Teledyne Isco (U.S.A.) - MPLC, HPLC, Fraction Collector, Syringe Pump
6. Zeton (U.S.A.) - Lab Scale Systems, Pilot Plants
7. SSI HPLC (U.S.A.) - HPLC Pump
8. Eldex Laboratories (U.S.A.) - MicroPro Pumps, ReciPro Pumps, Accessories
9. Chemyx (U.S.A.) - Syringe Pumps
10. Glas-Col (U.S.A.) - Heating mantle, Tape, Shaker, Stirrer, Rotators
11. VWR International (U.K.) - Water Purification Systems, General Lab
12. Amar Equipments (India) - Autoclaves, High Pressure-Temperature Stirred, Chemical Reactor
13. KGW Isotherm (Germany) - Dewar Flasks, Cold Traps, Tempering Vessels
14. Cecil (U.K.) - UV/Visible Spectrophotometers
15. Ultrawave (U.K.) - Cleaning Baths, Tank
16. MRC (Israel) - General Lab Equipment
17. DLab - Centrifuges, Pipettes, Thermo Control, Mixing
18. Optimize Technology (U.S.A.) - HPLC spare part : seal, piston, check valves
19. VICI Jour research - HPLC Fitting, Tubing, Valve, Check Valve
20. IDEX Health Science (U.S.A.) - HPLC Fitting, Tubing, Check Valve, Union

บริษัท ฟอ์จูน ไชแอนทิฟิค จำกัด

7/25 ถนนพหลโยธิน แขวงอนุสาวรีย์

เขตบางเขน กรุงเทพฯ 10220

Tel. 02-9861250-4 Fax. 02-9861255

Email : sales@fortunesci.com

Website : www.fortunesci.com

บริษัท ฟอ์จูน ไชแอนทิฟิค จำกัด เป็นตัวแทนจำหน่ายเครื่องมือวิเคราะห์ทางวิทยาศาสตร์คุณภาพสูง
ได้รับการยอมรับให้ใช้งานในหน่วยงานต่างๆ ทั้งในและต่างประเทศ ทั้งภาครัฐและเอกชน
ท่านสามารถมั่นใจได้ ด้วยประสบการณ์ทำงานยาวนานกว่า 20 ปี
พร้อมให้บริการด้วยทีมงาน และนักเคมีที่ผ่านการฝึกอบรมจากต่างประเทศ

โทรหาเราก่อนตัดสินใจซื้อ

แล้วท่านจะได้รับสินค้าและบริการในราคาที่พึงพอใจ